OVENS and INCUBATOR

User Manual





Model	Description	Temperature range
TCN-50 Super	Natural convection oven 50L (maximum volume)	Room temperature from + 5 °C to + 300 °C
TCN-115 Super	Natural convection oven 115L (maximum volume)	Room temperature from + 5 °C to + 300 °C
TCN-200 Super	Natural convection oven 200L (maximum volume)	Room temperature from + 5 °C to + 300 °C
TCF-50 Super	Air forced oven 50L (maximum volume)	Room temperature from + 10 °C to + 300 °C
TCF-120 Super	Air forced oven 120L (maximum volume)	Room temperature from + 10 °C to + 300 °C
TCF-200 Super	Air forced oven 200L (maximum volume)	Room temperature from + 10 °C to + 300 °C
TCF-400 Super	Air forced oven 400L (maximum volume)	Room temperature from + 10 °C to + 300 °C

Forced air and natural convection **ovens**, multifunctional with microprocessor temperature controller.

Forced air and natural convection **incubator**s, multifunctional with microprocessor temperature controller.

Model	Description Temperature range	
		Room temperature from + 5 °C to + 70 °C
ICN-35 Super	Natural convection incubator 35L (maximum volume)	Room temperature from + 5 °C to + 70 °C
ICN-55 Super	Natural convection incubator 55L (maximum volume)	Room temperature from + 5 °C to + 70 °C
ICN-120 Super	Natural convection incubator 120L (maximum volume)	Room temperature from + 5 °C to + 70 °C
ICN-200 Super	Natural convection incubator 200L (maximum volume)	Room temperature from + 5 °C to + 70 °C
ICF-55 Super	Air forced incubator 55L (maximum volume)	Room temperature from + 5 °C to + 80 °C
		(sterilisation special program at 130°C)
ICF-120 Super	Air forced incubator 120L (maximum volume)	Room temperature from + 5 °C to + 80 °C
		(sterilisation special program at 130°C)
ICF-200 Super Air forced incubator 200L (maximum volume) Room temperature from + 5 °		Room temperature from + 5 °C to + 80 °C
		(sterilisation special program at 130°C)
ICF-400 Super	Air forced incubator 400L (maximum volume)	Room temperature from + 5 °C to + 80 °C
		(sterilisation special program at 130°C)

Producer:

Drafting by the agent:

Sozhou Being Medical Device CO., LTD NO.108 Gongxiang RD Qiandeng Town Kunshan China

Giorgio Bormac s.r.l. Via della Meccanica, 25 41012 Carpi (MO) P. Iva 02309180368

Tel. +39 059 653274 Fax +39 059 653282 Email <u>info@giorgiobormac.com</u>

Index

1.	SAFETY INFORMATION
•	Definitions of warning words and symbols4
•	Reporting terms:
•	Warning symbols:
•	Pictograms
•	Danger symbols5
•	Prohibition signs
•	Symbols of obligation
2.	GENERAL SAFETY INSTRUCTIONS
3.	CE MARKING DATA
4.	CONTENT OF PACKAGE9
5.	TRANSPORTATION
•	Instructions for safe transportation10
•	Transportation of an already used oven/incubator10
6.	CONSERVATION
7.	FIRST INSTALLATION
•	Getting started10
8.	INSTRUMENT PARTS11
•	Display and commands12
9.	TECHNICAL SPECIFICATIONS
10.	OPERATING MODE14
•	Natural convection ovens / incubators14
•	Forced air ovens14

•	Forced air incubators15
11.	OPERATION15
•	Switching on the Instrument15
•	Basic Mode (PROG 0) or with Programs (PROG 18)15
•	Program recall and choice
•	Programming17
•	Start/stop of a program18
12.	USER SETTING MENU
•	Date and time setting20
•	Data recording frequency20
13.	INTRODUCTIONS OF SAMPLES INTO OVEN/INCUBATOR
•	Sample loading21
•	Temperature limit for sample protection21
•	Restart mode after absence of power energy22
•	Temperature limit for over-temperature alarm
•	Temperature safety device
•	Temperature offset on single point, on entire range, on room temperature sensor22
•	Delayed start function23
•	Repeat of the selected program23
14.	DOWNLOAD DATA AND USB-KEY23
15.	CLEAN AND MAINTENANCE24
16.	WARRANTY25
17.	DISPOSAL OF ELECTRONIC EQUIPMENT

1. Safety information

• Definitions of warning words and symbols

This manual contains extremely important safety information, to avoid personal injury, damage to the instrument, malfunctions or incorrect results due to failure to comply with them. Read entirely and carefully this manual and be sure to familiarize with the tool before starting to work with it. This manual must be kept near to the instrument, so that the operator can consult it easily, if necessary. Safety provisions are indicated with warning terms or symbols.

• Reporting terms:

ATTENTION/WARNING/DANGER	for a medium-risk hazardous situation, which could lead to serious injury or death, if not avoided.
WARNING	for important information about the product.
NOTE	for useful information about the product.

• Warning symbols:



DANGER

This symbol indicates an **imminently hazardous** situation, which, if not avoided, could result in death or serious (irreversible) injury.



WARNING

This symbol indicates a potential hazardous situation, which, if not avoided, could result in death or serious (irreversible) injury.



ATTENTION

This symbol indicates a potential hazardous situation, which, if not avoided, could result in medium or minor injuries (reversible).



ADVICE

This symbol draws attention to possible damage to the instrument or instrumental parts.



NOTE

This symbol highlights further information and tips.

• Pictograms

In this manual, there are various symbols identifying dangers, prohibitions, and obligations as illustrated below.

• Danger symbols



• Prohibition signs

Do not wet with water	
	Do not wet with water

Symbols of obligation

Disconnect the instrument from the power supply by pulling the plug	
Eye protection must be used	

2. General safety instructions

If the oven/incubator is not installed, operated, cleaned, adjusted or set up correctly, there is a risk of malfunction which could lead to personal injuries and material damages to the instrument and samples. Therefore, the oven/incubator must only be installed operated, cleaned, adjusted and set up by qualified personnel.



	WARNING
	Fire hazard The ovens/incubators must not be used if the inspection of thermostat with safety class 2 has failed.
 If the safety thermostat test fails, immediately stop using the oven/incubator, remove to plug from the power supply and contact your dealer for the necessary repairs. Always place the instrument on a working surface that is resistant to a temperature of 1 °C. 	
	igtriangle Do not put anything under the instrument (paper, plastic film, etc.).
	Always connect the instrument only to a power supply with a fuse of at least 10A. Follow the recommendations of your local power supply company.
	 Danger of burns The air intake cover on the back of the instrument becomes hot and must not be touched during operation of the oven.
	 Risk of injury and breakage Always place the instrument only on surfaces that can support its weight.
	Overturning hazard and risk of injury
	○ Never stack more than 2 ovens/incubators on top of each other.
	Always fix 2 stacked ovens with the fixing plates provided.
	Risk of injury, Risk of slipping or overturning and risk of damage to the instrument
	 The instrument must be lifted by 2 people. The instrument must be transported in its original packaging only.
	 The instrument must be transported in its original packaging only. The instrument must always be lifted from below with mechanical tools (e.g., forklift truck) together with the supporting pallet.
	 The instrument must not be lifted directly from below with mechanical tools without supporting pallets (e.g., forklift truck). The instrument must not be lifted or dragged by pulling the door.

3. CE Marking Data

ArgoLab instruments are manufactured in compliance with Directive 2006/42/EC and with the relevant Community Directives applicable at the time of its placing on the market (fac-simile below).

and Annex VI - Directive 2011/65/UE (RoHS)	SUZHOU BEING MEDICAL DEVICE.CO.,LTD	DECLARATION OF CONFORMITY UE In accordance with Annex II A - Directive 2006/42/CE Annex IV - EMC Directive and Annex VI - Directive 2011/65/UE (RoHS)	CE
--	-------------------------------------	--	----

No. ISETC.002420200624

Manufacturer's Name	: SUZHOU BEING MEDICAL DEVICE CO., LTD
Manufacturer's Address	: NO. 108 GONGXIANG RD QIANDENG TOWN, KUNSHAN CHINA
	Tel: +86-21-56633709
	Email: JILL.SHEN@BLUEPARD.COM
Object of Declaration:	: FORCED AIR INCUBATORS

This declaration of conformity is issued under the sole responsibility of the manufacturer.

Product names:	
Product description	FORCED AIR INCUBATORS
Model:	BI-120FL, BI-120F, BI-200FL, BI-200F, BI-400FL, BI-400F
Serial Number:	from s/n xxxxxxxxx to xxxxxxxxxx
Product options:	This declaration covers all options of the above products

The object of the declaration describe above complies with the essential requirements of the . following applicable European Directives, and carries the CE marking accordingly:

EMC directive: 2014/30/UE	Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility.
RoHS Directive 2011/65/EU	Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.
LVD Directive: 2014/35/UE	Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to the making available on the on the market of electrical equipment designed for use within certain voltage limits Text with EEA relevance.
Machinery Directive 2006/42/EC	DIRECTIVE 2006/42/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 17 May 2006 on machinery, and amending Directive 95/16/EC (recast)

• and conforms with the following standards:

EN 61010-1:2010+A1:2019

EN 61326-1:2013

EN 61000-3-2:2014

EN 61000-3-3:2013

EN 60204:2018

EN ISO 12100:2010

NAME AND ADDRESS OF THE PERSON AUTHORISED TO COMPILE THE TECHNICAL FILE

Giorgio Bormac S.r.l. - Via della Meccanica, 25 41012 Carpi (MO) - ITALY

Signed for and on behalf of	name, surname	
Place	gg/mm/aaaa	
SHANGHAI	SIGNATURE	

Fac-simile of the CE marking plate:

A	ARGO LAB	Made in P.R.C.
er	S/N 201148965	Date: 2020. 11
ce	Temp.Range RT	+ 5 °° ~ 70 °°
CE	Volts 220V/50Hz	Watts 600W
kat TOE Gongkang Rd. Kumithen Chine	Model BIT-200/	/ICN-200 Plus
being	Name In	cubator

4. Content of package

This instrument is delivered complete with the following parts:

- 1. n. 2 stainless steel wire shelves
- 2. n. 4 brackets for shelves
- 3. power supply cable
- 4. fuses
- 5. User manual
- 6. USB-Key

5. Transportation

• Instructions for safe transportation

 Risk of injury, of slipping or overturning and risk of damage to the instrument The instrument must be lifted by 2 people. The instrument must be transported in its original packaging only. The instrument must always be lifted from below with mechanical tools (e.g., forklift truck) together with the supporting pallet.
 The instrument must not be lifted directly from below with mechanical tools without supporting pallets (e.g., forklift truck). The instrument must not be lifted or dragged by pulling the door.

• Transportation of an already used oven/incubator



- Switch off the Argolab oven/incubator by pressing the main switch.
- Remove the power plug from the socket.
- Remove the shelves.
- Clean the Argolab oven/incubator and its shelves (see chapter 13 on page 31).
- Dry the inside of the Argolab oven/incubator and the shelves.
- Wrap the shelves with bubble wrap.
- Pack the shelves in their original packaging and place them in the Argolab oven/incubator.
- Pack the Argolab oven/incubator in its original packaging.
- Take care that the Argolab oven/incubator does not get wet during transport.
- During transport, maintain the permitted room temperature (from -10 °C to 60 °C).

6. Conservation

- Store the Argolab oven/incubator only in closed, dry rooms.
- The permitted storage temperature is -10 °C to 60 °C. The maximum permitted storage humidity is 85% RH without condensation.

7. First installation

• Getting started

The oven should be installed in following conditions:

- Dry, clean, stable worktable with a flat horizontal surface and heat resistant.
- At least 30 cm free around the instrument.
- Room temperature between 5 °C and 40 °C and relative humidity maximum of 85%.
- Power supply socket with earth connection.
- Power supply of 220/240 V 50 Hz.





8. Instrument parts



• Display and commands



Figure 3

The touch screen display is capacitive, so you can use it with each finger, wearing latex gloves, with pens with rubber pad, with generic pens or pencils.

NOTE: it is however recommended not to use the latter in order to avoid marking the display glass. The different colours of the display, icons and symbols have been used in the screens to make them easy to use and understandable:



COMMAND COLOUR	DEFINITION
White	Editable value or parameter
Grey	Not editable value or parameter
Green	START button, OK or ENTER button to confirm
Red	STOP button, indication of current temperature (> 20,0°C)
Orange	ESC button, operating parameters of the program steps (PROG 18)
Yellow	"Alarm" icon
Blue	indication of current temperature (\leq 20,0°C)

9. Technical specifications

Natural convection ovens	TCN-50 Super	TCN-115 Super	TCN-200 Super
Volume	50 l	115	200
Max temperature/ Resolution	+300/0,1°C	+300/0,1°C	+300/0,1°C
Homogeneity at 150°C	± 3,5°C	± 3,5°C	± 4,0 °C
Temperature variation at 150°C	± 0,5°C	± 0,5°C	± 0,7°C
Heating time at 150°C	16 min.	18 min.	20 min.
Timer	99:59 hh:min e ∞	99:59 hh:min e ∞	99:59 hh:min e ∞
Safety class	3.1	3.1	3.1
Power supply/power	230 V / 1000 W	230 V / 1900 W	230 V / 2100 W
Internal dimensions (L*A*P)	400 x 420 x 330 mm	520 x 495 x 450 mm	650 x 640 x 495 mm
Shelves number (standard/max)	2/5	3/6	3/9
Distance between shelves	50 mm	50 mm	50 mm
Max loading of shelves	15 Kg	20 Kg	20 Kg
External dimensions (L*A*P)	690 x 635 x 470 mm	815 x 750 x 600 mm	940 x 905 x 660 mm
Weight	53 Kg	74 Kg	103 Kg

Forced air ovens	TCF-50 Super	TCF-120 Super	TCF-200 Super	TCF-400 Super
Volume	50	1201	200	400
Max temperature/ Resolution	+300/0,1°C	+300/0,1°C	+300/0,1°C	+300/0,1°C
Homogeneity at 150°C	± 2 %	± 2 %	± 2 %	± 2 %
Temperature variation at 150°C	± 0,3°C	± 0,3°C	±0,4°C	± 0,5°C
Heating time at 150°C	20 min.	24 min.	30 min.	50 min.
Timer	99:59 hh:min e ∞			
Safety class	3.1	3.1	3.1	3.1
Power supply/power	230 V / 980 W	230 V / 1900 W	230 V / 2400 W	230 V / 3200 W
Internal dimensions (L*A*P)	400 x 415 x 310 mm	520 x 530 x 435 mm	645 x 650 x 495 mm	1000 x 800 x 500 mm
Shelves number (standard/max)	2/5	3/7	3/9	3/10
Distance between shelves	50 mm	50 mm	50 mm	50 mm
Max loading of shelves	15 Kg	20 Kg	20 Kg	20 Kg
External dimensions (L*A*P)	690 x 635 x 570 mm	810 x 750 x 690 mm	945 x 870 x 755 mm	1285 x 1060 x 750 mn
Weight	54 Kg	74 Kg	103 Kg	160 Kg

Natural convection ICN-16 Super		ICN-35 Super	ICN-55 Super	ICN-120 Super	ICN-200 Super
incubators					
Volume	16	35 I	55 l	1201	200
Max temperature/	+70/0,1°C	+70/0,1°C	+70/0,1°C	+70/0,1°C	+70/0,1°C
Resolution					
Homogeneity at 37°C	± 0,4 °C	± 0,4 °C	± 0,5 °C	± 0,5 °C	± 0,5 °C
Temperature variation at 37°C	± 0,3°C				
Heating time at 37°C	18 min.	22 min.	25 min.	30 min.	35 min.
Timer	99:59 hh:min e ∞				
Safety class	2	2	2	2	2
Power supply/power	230 V / 85 W	230 V / 125 W	230 V / 250 W	230 V / 350 W	230 V / 600 W
Internal dimensions (L*A*P)	270 x 230 x 255 mm	360 x 300 x 320 mm	400 x 360 x 385 mm	520 x 460 x 500 mm	610 x 600 x 575 mm
Shelves number (standard/max)	2/3	2/5	2/5	3/7	3/9
Distance between shelves	25 mm	30 mm	50 mm	50 mm	50 mm
Max loading of shelves	5 Kg	7,5 Kg	10 Kg	10 Kg	10 Kg
External dimensions (L*A*P)	530 x 370 x 400 mm	620 x 440 x 460 mm	660 x 500 x 545 mm	780 x 610 x 645 mm	875 x 755 x 710 mm
Weight	23 Kg	33 Kg	42 Kg	61 Kg	77 Kg

Forced air incubators	ICF-55 Super	ICF-120 Super	ICF-200 Super	ICF-400 Super
Temperature variation at 37°C	57	120	200	400
Heating time at 37°C	+80/0,1°C	+80/0,1°C	+80/0,1°C	+80/0,1°C
Timer	± 0,3 °C	± 0,4 °C	± 0,4 °C	± 0,5 °C
Safety class	± 0,1°C	± 0,1°C	± 0,2°C	± 0,3°C
Power supply/power	30 min.	40 min.	45 min.	55 min.
Internal dimensions (L*A*P)	99:59 hh:min e ∞			
Shelves number (standard/max)	3.1	3.1	3.1	3.1
Distance between shelves	230 V / 350 W	230 V / 600 W	230 V / 700 W	230 V / 1500 W
Max loading of shelves	400 x 415 x 350 mm	520 x 530 x 435 mm	645 x 650 x 495 mm	1000 x 800 x 500 mm
External dimensions (L*A*P)	2/5	2/7	3/9	2/10
Weight	50 mm	50 mm	50 mm	50 mm
Temperature variation at 37°C	20 Kg	20 Kg	20 Kg	20 Kg
Heating time at 37°C	690 x 650 x 620 mm	810 x 750 x 690 mm	945 x 870 x 755 mm	1285 x 1060 x 750 mm
Timer	56 Kg	74 Kg	103 Kg	160 Kg

10. Operating mode

• Natural convection ovens / incubators

The instruments of the TCN and ICN series are with natural convection.

This means that, in the internal chamber of the instrument, the heat is propagated through the natural convective motions created by the heat exchange between cold and hot air. In ArgoLab natural convection instruments, there are special manual valves aimed at the correct recirculation of air inside the chambers of ovens and incubators. **NOTE**: ArgoLab instruments are supplied with the valves open, it is advisable not to close them in order not to affect the performance of the instrument.

NOTE: depending on the model, the lower valves will be present or not.



• Forced air ovens

The instruments of the TCF series are with forced ventilation.

This means that, in the internal chamber of the instrument, the heat is distributed homogeneously through the designated fan. In ArgoLab forced ventilation ovens (TCF series), there is an adjustable manual valve (dedicated to incoming cold air) aimed at the exchange of air inside the chamber.

NOTE: ArgoLab ovens are supplied with open valve, it is recommended not to close it completely in order not to affect the performance of the instrument.

NOTE: in the TCF 400 model there are 2 valves for the discharge of hot air (located at the top) and 2 valves for the loading of cold air (located at the bottom), for each discharge valve there is a fan.





• Forced air incubators

The instruments of ICF series are with forced ventilation. This means that, in the internal chamber of the instrument, the heat is distributed homogeneously through the designated fan.



11. Operation

• Switching on the Instrument

Connect the power cable to a grounded socket. Switch on the instrument using the ON/OFF button (see Figure 1). The button and the display light up. The display shows the sign "Super" in the first screen and the installed software version.

• Basic Mode (PROG 0) or with Programs (PROG 1...8)

The instrument can work in 2 modes:

- Basic (PROG 0) single step operation program.
- Programs (PROG 1...8) 8 recordable programs, each of which consists of 8 steps.

In both versions, it is always possible to set the temperature, the timer and the fan speed (where present) for each single work step.

Depending on the program you are in, the display will show one of the following main screens: **PROG 0**

PROG 1...8







Figure 7 – Main screen PROG 0

• Program recall and choice

To choose the program you want to use, press the key highlighted in Figure 7 or 8. It is in the same position, as mentioned, what changes is only the content depending on the program you are using. One of the 2 screens (Figure 9 or 10) will then appear.



ESC)	сно	OSE	PRO	GRAN	8 1	(ОК
Step	1	2	3	4	5	6	7	8
-	300. 0	100.0	150.0	200.0	250. 0	300.0	30.0	30. 0
	5999	120	150	180	0	120	0	0
20	100%	100%	100%	75%	100%	100%	75%	75%

Figure 9

Figure 10

As mentioned, the PROG has a single step, hence the indication "FIX MODE–NO PROGRAM AVAILABLE" (Figure 7), while programs from one to eight have eight steps each (Figure 8).

NOTE: in the screen of Figure 10, the operating parameters of the various steps of the selected program are visible in read-only mode (grey colour), but not editable. To be able to modify them, it is necessary to follow the instructions in the paragraph "Modifying operating parameters PROG 1...8".



Using the key highlighted in yellow, it possible to access the numeric keypad in Figure 11 and select the desired program.





Press ENTER to confirm the entered value. The display return to the previous screen where it is necessary to press OK to validate and then return to the main screen.

NOTE: it is possible to return to the previous screen without saving any changes by pressing ESC.



Programming •

• Modification of PROG 0 operating parameters:

From the main screen of PROG 0 mode (Figure 12), press the highlighted keys to modify respectively from left to right, the temperature, the fan speed (if present) and the timer. The screens in Figure 13, Figure 14, Figure 15 will therefore appear.





Figure 14

To increase or decrease the temperature and time values, use the + und - keys, while for the fan (if present), it is possible to choose directly between three speeds: High (100%), Medium (75%) e Low (50%). In all the screen, it is necessary to confirm the value set with the OK button.

NOTE: it is possible to return to previous screen without saving changes by pressing ESC.

• Modification of operating parameters PROG 1...8:

From the main screen (Figure 16), press the PROG key to access the menu in Figure 17.





In the screen of Figure 17, select the program to be modified by pressing the key highlighted in yellow and setting the corresponding digit using the numeric keypad. Confirm with ENTER button.

For each single program step (from 1 to 8), it is possible to set the temperature (degrees centigrade), the timer (minutes) and the speed of the fan (if present), by clicking respectively on the keys highlighted in green, red and pink of the same screen.

The respective numeric keypads will appear for the temperature and the timer, marked by the thermometer icon (Figure 18) and the clock (Figure 19), indicating the program and the step being edited. Once desired value has been set, confirm with ENTER or return to the previous screen without saving the changes by pressing ESC.



thout saving the changes by pressing Ese.
SET PROGRAM 8 STEP 1
RANGE [0 ~ 5999]
67890 -

Figure 18

Figure 19

To change the fan speed (if present), on the other hand, simply click on the corresponding button (pink in Figure 17) and the value will change in sequence between 100% = High, 75% = Medium and 50% = Low.

REPEAT THE PREVIOUS INSTRUCTIONS FOR EACH STEP YOU WANT TO PROGRAM

NOTE: If you do not want to use all the STEPS of the program you are storing, you must force the instrument to terminate the program itself. To do this, simply set a time equal to zero in the step following the last one you want to use.



EXAMPLE of Figure 17: If the last work step to be used is the third, simply set the timer equal to 0 in the fourth step, thus requiring the instrument to stop after the third step.

• Start/stop of a program

To start or stop a program, press respectively the START (Figure 20) or STOP (Figure 21) button.



Figure 20

Figure 21

The instrument therefore begins the set work cycle which can consist of one or more steps. **NOTE**: the timer countdown starts only when the set temperature is reached. For accuracy when T real= T set \pm 0,3°C.



Looking at Figure 20, for example, the countdown of 99 hours and 59 minutes will start when the instrument reaches 300 ± 0,3° °C. At any time, it is possible to check the work step and the number of the program in which you are (in green Figure 22 and Figure 23), as well as the possible delay (Delay) that has been set and the number of repetitions of the performed program (Cycle), highlighted in blue and in Figure 22 and 23.



Figure 22

Figure 23

NOTE: it is always possible to consult the various parameters set in the steps of the program in use by clicking on the green button of Figure 22 and Figure 23 but be careful to press ESC only and exclusively to return to the main screen. If you press the OK key, the instrument will interpret the action as confirming the choice of a new program, therefore it will forcefully terminate the work cycle in progress. When the program ends or is stopped manually, the word "End" appears on the main screen (Figure 24) together with an intermittent acoustic signal, which ca be silenced by pressing the word "End" itself.



Figure 24

User setting menu 12.

By pressing the gears icon in the main display (Figure 25) and then the USER button (Figure 26), you access the user settings menu screens (Figure 27 and Figure 28).

⇔ memory	01 Oct 2018	10:33			
23.7 ℃	00:00) hh:mm	USER	SERVICE	ADMIN
SET SET 0.0 C		STEP 1 PROG 1		L	
START End	PROG	\$ °	ESC		*
Figure 25					Figure 26

Figure 25



Figure 27

Figure 28

Date and time setting

By clicking on any button of date or time, you access the TIME menu (Figure 29).



Figure 29

Set the date and time using the numeric keypad, bearing in mind that the format is as follows:

20yy-mm-dd hh:mm:ss

NOTE: It is necessary to set all parameters (date and time) each time.



Data recording frequency

Click on the button highlighted in Figure 30 to modify the time frequency with which the instrument records the data of work cycle.





Using the numeric keypad, set a time (in minutes) from 1 to 255. Press ENTER to confirm or ESC to return to the previous screen without saving.

NOTE: the memory can hold up to 2000 records. Once this limit has been reached, the instrument will begin to overwrite the memory itself with new data, thus erasing the previous ones starting from the oldest. The recording frequency therefore also defines the recording capacity in terms of duration (eg frequency = 60 minutes x 2000 data x 60 minutes = 2000 hours = 83 days).

NOTE: if a USB key is connected to the instrument, the recordings also take place directly on it as well as in the instrument's memory. The memory limit of 2000 data remains, but the key itself will act as secondary memory and therefore, it is possible to carry out a continuous recording well over 2000 data.



13. Introductions of samples into oven/incubator

Explosion hazard and risk of death
○ Never introduce materials into the instrument that are explosive or flammable at the selected operating temperature.
⊘ Never introduce into the instrument's materials containing flammable or explosive solvents.
\odot Never introduce into the instrument materials that by sublimation or pyrolysis create
flammable materials at the selected operating temperature.
Poisoning and death hazard
⊘ Never introduce materials into the instrument which could create poisonous gases.
\odot Never introduce materials into the instrument that can react with moisture and produce
explosive gases.

• Sample loading

To obtain optimal air circulation inside the chamber of Argolab ovens/incubators, it is recommended to leave empty spaces between the samples (see Figure 31).

For correct convection of the samples, it is recommended not to put them in contact with the walls of the ArgoLab oven/incubator.





• Temperature limit for sample protection

The instrument expects to be able to limit the maximum working temperature (Tm) for the protection of samples from an incorrect setting of the temperature of the heating cycle. The "Tm" parameter (max temperature) and the maximum value foreseen by that type of instrument appear on the display at the top left of Figure 28 (different between oven and incubator). Set the maximum temperature value that you do not want the instrument to exceed during operation by pressing on the box and using the numeric keypad. Press ENTER to confirm or ESC to return to the previous screen without saving.

Example: if the temperature set for the heating cycle is 100 °C and a limit temperature (Tm) of 70°C is set, the instrument will try to reach the temperature set during the parameter setting (100°C), even if it is higher than the limit set in this submenu (Tm). When 70 °C is reached, the instrument goes into alarm by emitting an intermittent acoustic signal (which can be silenced by pressing the Alarm button) and the heating element is no longer powered until the temperature drops below that limit.

NOTE: to determine the correct "**tm**" value, the natural and inevitable initial temperature peak that the Argolab oven/incubator will have during thermostatting must be considered.



Application example: If the temperature set for the heating cycle is 100 °C and a limit temperature (tm) of 70 °C is set, the instrument will attempt to reach the temperature indicated during parameter setting (100 °C), even if it is higher than the limit temperature set in this submenu (tm). When 70 °C is reached the instrument goes into alarm with an intermittent acoustic signal (can be silenced by pressing any key) and the heating element is no longer powered until the temperature falls below the limit temperature ("tm").

NOTE: the instrument will always attempt to reach the temperature set for the heating cycle and consequently, as long as it is higher than the limit temperature, the device will go into an over-temperature alarm.



• Restart mode after absence of power energy

The mode in which the instrument starts working again after a power failure (Po) can be set:

VALUE Po	DESCRIPTION
0	When the power supply returns, the instrument does not automatically resume the heating cycle but must be restarted manually.
1	When the power supply returns, the instrument automatically resumes operation from the start of the interrupted heating cycle.
2	When the power supply returns, the instrument automatically resumes operation from the precise point when the heating cycle was interrupted.

Set the desired value by pressing on the corresponding box (figure 28) and using the numeric keypad. Press ENTER to confirm or ESC to return to the previous screen without saving.

• Temperature limit for over-temperature alarm

It is possible for the user to set the temperature value beyond which the instrument goes into overtemperature alarm (AI).

NOTE: although it can be modified by the operator, this value is already set at the factory and is specifically calibrated to the type of instrument in question, natural/forced oven or natural/forced incubator. It is therefore advisable not to modify this value unless strictly necessary, as temperature



fluctuations above or below the set value, especially in natural convection models, are completely normal and therefore reducing the AL value excessively would risk causing the instrument to go into alarm frequently and unnecessarily. For any changes, click on the corresponding box (Figure 28) and use the numeric keypad. Press ENTER to confirm or ESC to return to the previous screen without saving.

• Temperature safety device

Each ArgoLab oven/incubator instrument has an electronic overtemperature limiter (Protection class 2 according to the technical standard DIN 12880).

ArgoLab TCN and TCF series ovens and ICF series forced ventilation incubators have an additional adjustable protection (electromechanical type).

The electromechanical safety device (Class 3.1 according to the technical standard DIN 12880) is installed inside the left side panel of the instrument.

• Temperature offset on single point, on entire range, on room temperature sensor

The instrument provides the possibility for the user to set the offset values, i.e. calibration, on a temperature single point (**Pb**), on entire range (**Pk**) and on room temperature one (**Pa**).

NOTE: although they can be changed by the operator, these values are already factory set and perfectly calibrated with Accredia certified and traceable measuring instruments. It is therefore advisable not to modify these values unless strictly necessary, for example if, by checking with a digital and certified thermometer inconsistencies are detected between the reading of the temperature values of the instrument.



thermometer, inconsistencies are detected between the reading of the temperature values of the instrument and those detected by the thermometer itself.

PARAMETER	DESCRIPTION
Pb	By modifying this parameter, it is possible to correct the reading of the PT100 temperature sensor inside the instrument to a single temperature point. The correction will therefore be referable to only one specific point.
Pk	By modifying this parameter, it is possible to correct the reading of the PT100 temperature sensor inside the instrument over the entire temperature range, i.e., it is possible to vary the inclination of the reading range of the sensor itself.
Ра	By modifying this parameter, it is possible to correct the reading of the PT100 room temperature sensor installed on the instrument (refrigerated versions only) to a single temperature point. The correction will therefore be referable to only one specific point.

To make any changes, click on the corresponding boxes (Figure 25) and use the numeric keypad. Press ENTER to confirm or ESC to return the previous screen without saving.

Delayed start function

The instrument provides for the possibility of setting a delay at the start of the operating cycle (Delay) from 1 to 9999 minutes. Set the desired value by pressing on the corresponding box (Figure 25) and using the numeric keypad. Press ENTER to confirm or ESC to return to the previous screen without saving. Once the delay has been set, by pressing the START key, the instrument starts the program, but does not immediately start heating. After the sett delay time has elapsed, the instrument starts heating and the set timer appears on the display.

Repeat of the selected program ۲

The instrument allows you to repeat the selected program from 1 to 99 times, function Cycle. **NOTE**: it is possible to set the continuous repetition of a program too by putting it in continuous "loop" by setting the parameter Cycle= 0. Set the desired value by pressing the corresponding box (Figure 25) and using the numeric keypad. Press ENTER to confirm or ESC to return to the previous screen without saving.

Download data and USB-Key 14.

The instrument is equipped as standard with a USB port for connecting a pen drive (Figure 32). Connecting a USB pen drive when the instrument is switched on (either with the work cycle active or inactive), a fully automatic download (without having to press anything) of all the data that the machine has in memory takes place.

NOTE: when the data download takes place, the internal memory is completely emptied and the data transferred to the USB pen drive. Therefore, in this way the internal memory is formatted and all data in it are deleted.

NOTE: in case a USB stick is left connected to the instrument, the recordings are made directly on it as well as in the instrument memory. The memory limit of 2000 data remains, but the stick itself will act as a secondary memory and it is therefore possible to make a continuous recording well beyond 2000 data.

The format of the files that are downloaded is .txt. The data are therefore completely "open", modifiable and transferable to normal computer applications, and can be processed independently by the operator, without the need for any dedicated software. The files are organized in folders divided by year (Y) and month (M) (Figure 33) and in turn within the various folders are ordered in days (D01, D02...). If more than one work cycle is executed on the same day, they will be named with the same date plus the progressive cycle number. In the example of (Figure 34) on day 3 four work cycles D03-0, D03-1, D03-2, D03-3 were run.



Figure 32







Unità USB (E:) > Y2018M09			NUM	DATE	TIME	REALTEMP	SET TEMP	
		~		1	2018-10-09	08:44	100.9	0.0
^	Nome			2	2018-10-09	08:45	88.2	0.0
	D03-0			3	2018-10-09	08:46	84.8	0.0
	D03-1			4	2018-10-09	08:47	74.9	0.0
	D03-2			5	2018-10-09	08:48	63.8	0.0
	D03-3			6	2018-10-09	08:49	55.6	0.0
	D06-0			7	2018-10-09	08:50	48.0	0.0
	D15-0			8	2018-10-09	08:51	42.5	0.0
	D16-0			9	2018-10-09	08:52	39.5	0.0
	D17-0			10	2018-10-09	08:53	39.2	40.0
	D25-0			11	2018-10-09	08:54	38.3	40.0
	D26-0			12	2018-10-09	08:55	39.7	40.0
	D26-1			13	2018-10-09	08:56	40.7	40.0
	D27-0			14	2018-10-09	08:57	41.0	40.0
	D27-1			15	2018-10-09	08:58	41.0	40.0
	D28-0			16	2018-10-09	08:59	40.7	40.0
	020-0			17	2018-10-09	09:00	40.7	40.0
Figu	re 34		Figure 35	18	2018-10-09	09:01	40.7	40.0
				19	2018-10-09	09:02	40.6	40.0
				20	2018-10-09	09:03	40.3	40.0

As mentioned, the .txt files are completely open and contain the main working data of the machine: registration number (NUM), date (DATE), time (TIME), detected temperature (REAL TEMP) and set temperature (SET TEMP), see Figure 35. They are already organized in tabular form and can be transferred to other computer applications in the form of text or table of values.

15. Clean and maintenance

A correct maintenance and cleaning of the instrument ensure that it remains in good condition.

The internal chamber of the instrument is made of stainless steel, so it can be cleaned with any detergent provided, that is not aggressive and/or corrosive.



It is recommended that internal and external surfaces are cleaned with a normal all-purpose cleaner sprayed onto a soft dampened cloth, so that it is not used in concentrated form. Before proceeding with cleaning or decontamination, the user must ensure that the method used does not damage the instrument.

ATTENTION							
	Danger of corrosion – Damage to the instrument						
	arnothing DO NOT use cleaning agents containing acids or halides.						
	arnothing DO NOT use neutral cleaning agents on other surfaces						
	(e.g., on galvanised parts of the hinges or the rear wall of the housing).						
	Eye contact – Eye damage caused by chemical burns						
	arnothing DO NOT discharge into the sewage system.						
	Wear safety goggles.						

IMPORTANT:

If the instrument is to be sent for service, it should be properly cleaned and possibly decontaminated from pathogens. It is also advisable to return the instrument in its original packaging to the repair service and if this is not possible to pack it adequately for transport. Any damage caused by incorrect shipment will not be covered by warranty.

16. Warranty

Under normal use this instrument is guaranteed for a period of 24 months from the date of purchase.

The warranty is valid only if the product purchased remains original. It does not apply to any product or parts thereof that have been damaged due to incorrect installation, improper connection, misuse, accident or abnormal operating conditions. No liability is accepted for damages caused by improper use, lack of maintenance and unauthorised modifications.

17. Disposal of electronic equipment



This equipment is subject to regulations for electronic devices. Dispose of in accordance with local regulations.